



Intelligent
health systems

TALKING

Forehead Thermometer

PROFESSIONAL ACCURACY



Make an intelligent choice...

WHY CHOOSE US?

Our thermometer takes the guesswork out of accuracy.

No messy probes, or invasive procedures.

Just place the thermometer in the centre of the forehead, press the measure button, and in 2 seconds you have an accurate reading, without any fuss.

Our thermometer will allow you to take a reading in the dark with its talking feature. It will also record the last 30 readings to see if there is any fluctuations for you to act upon. This ensures that you never have to wake up a child or administer medicine unnecessarily.

We believe that it is the best, and so will you.

WHY IS IT MORE ACCURATE THAN AN EAR THERMOMETER?

The major reason ear thermometers are considered inaccurate by medical professionals is due to the positioning of the probe in the ear canal being inconsistent, thus creating inconsistent readings and frequently missing fevers. The Forehead thermometer scans the forehead area for the temporal artery and it is almost impossible to miss the artery during a scan.



- Easy to use
- Fast and reliable
- Portable size
- Easy to read with large LCD screen
- Professional accuracy
- Temperature taking without pain or tears
- 12 months warranty

Available at all major department and independent stores, including leading pharmacies.

For more information, please email: info@intelligenthealthsystems.com.au or visit: www.intelligenthealthsystems.com.au